


Next Generation Multiplexing Platform

Luminex® xMAP®
INTELLIFLEX® System



The Life Science business
of Merck operates as
MilliporeSigma in the
U.S. and Canada.

Millipore®

Preparation, Separation,
Filtration & Monitoring Products

The Luminex® xMAP® INTELLIFLEX® system is a compact, flow-based, multiplex platform that simplifies your workflow by combining the reliable performance of xMAP® technology with advanced features to elevate your assay performance and empower assay development innovation.

By combining low- and high-plex capabilities, the xMAP® INTELLIFLEX® system gives

fast, reliable results.

Features of the xMAP® INTELLIFLEX® System:

- **Wide Dynamic Range**
Avoid repeat testing due to out-of-range values with the widest dynamic range of any xMAP® platform.
- **Compact Footprint**
Reduce lab bench requirements with a compact footprint that features an embedded PC and a touchscreen user interface.
- **User-Friendly and Flexible Software**
Quickly master the platform with easy-to-use software that allows you to optimize workflows, gather performance insights, and export data into common analysis tools.
- **Compatible with Existing MILLIPLEX® Assays**
Experience the versatility of xMAP® technology with an open architecture design, enabling you to choose from the most extensive menu of MILLIPLEX® assays and custom assay solutions.
- **Dual Reporter Functionality**
Now get more data from each target with the second reporter channel that allows you to save on time and precious sample.



Optimize Your Workflow

Quickly assemble plate configurations from saved components and start your run in as few as 20 seconds. Plate configurations can be set up with just a few clicks or can be imported from a CSV template.





Figure 1. xMAP® INTELLIFLEX® exemplary layout controls:

- 1 Define plate layouts, panels, or acquisition parameters in any sequence, and save modular components for use in different acquisitions.
- 2 Import plate layouts, sample lists, or panels from CSV.
- 3 Define replicates with intuitive copy-and-paste controls.
- 4 View run parameters at a glance.


Monitor Performance Fast

Quickly navigate between multiple data types for comprehensive, real-time run monitoring.


Bead Map



Doublet Discriminator Gate



MFI Heat Map



Well Status




Figure 2. xMAP® INTELLIFLEX® run monitoring data examples.

Easily Export Data for Analysis

Export results in one of two CSV formats—an xPONENT®-style format or an INTELLIFLEX® style data format—designed for flexibility and simplicity. The xPONENT®-style data format can be analyzed using our Belysa® Immunoassay Curve Fitting Software.





Figure 3. Export selection screen for xMAP® INTELLIFLEX® data format.

Get More Data with Dual Reporter Functionality

With the addition of a second reporter channel, you can acquire data about two parameters per target protein or nucleic acid, saving on time and precious sample volume. This is ideal for measuring multiple isotypes or configurations that target the same antigen.






Figure 4: Equivalent signal in both single and dual reporter modes. Phycoerythrin-labeled IgM detection antibody and streptavidin-Super Bright 436 + biotinylated IgG detection antibody were used to measure SARS-CoV-2 N antigen-specific IgM and IgG in single and dual reporter modes. When run together, no interference was observed between the dyes and MFI values were equivalent.

For Research Use Only.
Not For Use In Diagnostic Procedures.

Specifications

Performance	xMAP® INTELLIFLEX®	xMAP® INTELLIFLEX® DR-SE
Multiplexing Capability	Up to 500 individual analytes per well	Up to 500 individual analytes per well Up to 2 parameters per analyte
Typical Dynamic Range	≥ 5.5 logs	≥ 5.5 logs (RP1) ≥ 4.5 logs (RP2)
Plate Run Time		96-well plate in ~20 min 384-well plate in ~75 min
Sensitivity	Commercially available kits deliver sensitivity as low as sub-picogram level	
General		
Physical Dimensions	58.4 cm (23 in.) W x 61 cm (24 in.) D x 76.2 cm (30 in.) H	
Weight	54.4 kg (120 lbs)	
System Warmup	30 min	
Operating Temperature	15–30°C (59–86°F)	
Operating System	Microsoft® Windows® 10 IoT Enterprise LTSC	
Barcode Reader	For importing target values from the xMAP® INTELLIFLEX® Calibration and Performance Verification Kits	
Fluidics		
Sheath Pressure	8–13 psi (15 psi maximum)	
Sample Injection Rate	2 µL/sec	
Sample Uptake Volume	10–200 µL	
Optics		
Lasers	1 Classification laser 1 Reporter laser	1 Classification laser 2 Reporter lasers
Classification Laser Excitation	638 nm	
Reporter Laser 1 Excitation	532 nm	532 nm (RP1) 405 nm (RP2)
Detector	565–585 nm	565–585 nm (RP1) 421–441 nm (RP2)

DR-SE = dual reporter side eject

Ordering Information

Product Description	Cat. No.
xMAP® INTELLIFLEX® RUO System	INTELLIFLEX-RUO

Includes:

- xMAP® INTELLIFLEX® RUO Instrument
- xMAP® INTELLIFLEX® Software (Installed)
- xMAP® INTELLIFLEX® Quick Guide
- xMAP® INTELLIFLEX® Performance Verification Kit
- xMAP® INTELLIFLEX® Calibration Kit
- xMAP® Sheath Fluid PLUS, 20 L
- (3) Sample Probe Needles
- Sample Probe Height Adjustment Kit
- Power Cable, Country-Specific
- Off-Plate Reagent Block
- Keyboard and Mouse Combo
- Heater Block for 96-Well Microtiter Plates
- Heater Block for 384-Well Microtiter Plates
- Empty Cubitainer for Waste

xMAP® INTELLIFLEX® DR-SE	INTELLIFLEX-DRSE
---------------------------------	-------------------------

Includes:

- xMAP® INTELLIFLEX® DR-SE RUO Instrument
- xMAP® INTELLIFLEX® Software (Installed)
- xMAP® INTELLIFLEX® Quick Guide
- xMAP® INTELLIFLEX® Performance Verification Kit
- xMAP® INTELLIFLEX® Calibration Kit
- xMAP® Sheath Fluid PLUS, 20 L
- (3) Sample Probe Needles
- Sample Probe Height Adjustment Kit
- Power Cable, Country-Specific
- Off-Plate Reagent Block
- Keyboard and Mouse Combo
- Heater Block for 96-Well Microtiter Plates
- Heater Block for 384-Well Microtiter Plates
- Empty Cubitainer for Waste

Ancillary Products	
xMAP® INTELLIFLEX® Performance Verification Kit RUO, 20 Uses	IFX-PVER-K20
xMAP® INTELLIFLEX® Calibration Kit RUO, 20 Uses	IFX-CAL-K20
xMAP® Sheath Fluid PLUS RUO, 20 L	40-50021
xMAP® Sheath Concentrate PLUS, 1 L	40-50023
Belysa® Immunoassay Curve Fitting Software	40-122
BioTek® 405™ LS Magnetic 96-Well Washer	40-092AB

Product Description	Cat. No.
Instrument Services	
Platinum Service	IFLEX-PLATINUM
Gold Service	IFLEX-GOLD
Platinum Service (DR-SE model)	IFLEX-DRSE-PLAT
Gold Service (DR-SE model)	IFLEX-DRSE-GOLD
Validation Services	
IQ/OQ with Installation	IFLEX-IQOQ
Yearly OQ	IFLEX-OQ

To learn more about the xMAP® INTELLIFLEX® System, visit SigmaAldrich.com/intelliflex

Or, explore our MILLIPLEX® portfolio of immunoassays, at SigmaAldrich.com/milliplex



Millipore®

Preparation, Separation,
Filtration & Monitoring Products

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

MerckMillipore.com

To place an order or receive technical assistance in Europe, please call Customer Service:

France: 0825 045 645

Spain: 901 516 645 Option 1

Germany: 069 86798021

Switzerland: 0848 645 645

Italy: 848 845 645

United Kingdom: 0870 900 4645

For other countries across Europe, please call: +44 (0) 115 943 0840

Or visit: **[SigmaAldrich.com/offices](https://www.sigmaaldrich.com/offices)**

For Technical Service visit: **[SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice)**

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, Millipore, MILLIPLEX, and Belysa are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MK_BR8111EN Ver. 4.0
46477
02/2023